



PRODUCT FEATURES

SPIROMETRY FILTER TUBE "KIT" (white-orange)

For use during pulmonary lung function testing Spirometers: CareFusion - Jaeger

Description: **Filter:** barrier filtration for protection of patient, staff and equipment from bacterial and viral contamination.
Filter + Nose clip "CUSHION"

Code items:

- individually packaged:

30909-M1, Adult / Paediatric (30811 + 30851-M!)

Materials:

- housing:
- filter media:
- nose clip:
- sealed:

Polystyrene (HIPS) recyclable
Hydrophobic Polypropylene (TECHN)
PolyAcetato (POM) - Durecon M270
Ultrasonically

Kind of filtration:

Electrostatic / Mechanical

Flow resistance: (complies with ATS)

- Flow rate 12,00 L/sec

< 0,93 cm H₂O / (L/s)

Bacterial Filtration Efficiency:

99,9 %

Viral Filtration Efficiency:

99,9 %

Dead space of the filter:
(connector not included)

32,00 ml

Filtration area:

11,46 cm²

Exchange filtration area:

> 34,00 cm²

Weight of the filter:

18,00 g

Working procedure:

CE 0373
Disposable

Class IIa - Italian Superior Institute of Health

Quantity / box:

150 pcs.

TEST: NELSON Laboratories - USA -

NOTES: single patient use. Latex free.

Medical Devices Regulations - Directive 93/42/EEC

Any size of adaptors is available upon request.

MADA

Spirometry Filters

MADA S.r.l. Via Cava Trombetta, 5 - 20054 Segrate - Milan (Italy)

Tel. +39 (0)2.26924623 - Fax +39 (0)2.26924602

e-mail: info@mada-filters.it - website: www.mada-filters.it