



PRODUCT FEATURES

SPIROMETRY FILTER **TUBE-C “KIT”** (white-green)

For use during pulmonary lung function testing **Spirometers: Cosmed - Medisoft**

Description: **Filter:** barrier filtration for protection of patient, staff and equipment from bacterial and viral contamination.

Filter + Nose clip “CUSHION”

Code items:

- individually packaged: **30912-M1 Adult / Paediatric** (30814 + 30851-M1)

Materials:

- housing: Polystyrene (HIPS) recyclable
 - filter media: Hydrophobic Polypropylene (TECHN)
 - nose clip: PolyAcetato (POM) - Durecon M270
 - sealed: Ultrasonically

Kind of filtration:

Electrostatic / Mechanical

Flow resistance: (complies with ATS)

- Flow rate 12,00 L/sec < 0,93 cm H₂O / (L/s)

Bacterial Filtration Efficiency:

Flow rate 30 L/min 99,9 %

Viral Filtration Efficiency:

Flow rate 30 L/min 99,9 %

Dead space of the filter:
 (connector not included)

29,52 ml

Filtration area:

11,46 cm²

Exchange filtration area:

> 34,00 cm²

Weight of the filter:

18,00 g

Working procedure:

CE 0373 Class IIa - Italian Superior Institute of Health
 Disposable

Quantity / box:

150 pcs.

TEST: NELSON Laboratories - USA -

NOTES: single patient use. Latex free.

Medical Devices Regulations - Directive 93/42/EEC

Any size of adaptors is available upon request.

MADA

Spirometry Filters

MADA S.r.l. Via Cava Trombetta, 5 - 20054 Segrate – Milan (Italy)
 Tel. +39 (0)2.26924623 - Fax +39 (0)2.26924602
 e-mail: info@mada-filters.it - website: www.mada-filters.it